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序号	引物名称	引物	限制性内切酶	Nmoles	所加TE，目标浓度为100uM
1	Prime1_A_bio	5' Bio-CCTGACGAGCATCACAAA 3'			
2	Prime1_B_xbai	5' GCTCTAGA GCCATACCAAACGACGAG 3'	XbaI		
3	Prime2_A_saci	5' CGAGCTC CCTGACGAGCATCACAAA 3'	SacI		
4	Prime2_B_dig	5' dig-GCCATACCAAACGACGAG 3'			

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One Taq® 2X Master Mix with Standard Buffer

One Taq® 2X Master Mix with Standard Buffer is an optimized blend of Taq and Deep VentR™ DNA Polymerases ideally suited to routine PCR applications from a variety of templates including pure DNA solutions, bacterial colonies, and cDNA products. The 3'→5' exonuclease activity of Deep Vent DNA Polymerase increases the fidelity and robust amplification of Taq DNA Polymerase (1). The convenient master mix formulation contains dNTPs, MgCl2, buffer components and stabilizers, requiring only the addition of primers and DNA template for robust amplification.

Buffer Composition
20 mM Tris-HCl
22 mM KCl
22 mM NH4Cl
1.8 mM MgCl2
5% Glycerol
0.05% Tween® 20
0.06% IGEPAL® CA-630
0.2 mM dNTPs
25 units/ml OneTaq® DNA Polymerase
pH 8.9@25°C

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Binding capacity	10 μg
Elution volume	> 30 μl
Format	Tube
Fragment size	100 bp – 10 kb
Fragments removed	< 40mers
Processing	Manual
Sample type: applications	ssDNA or dsDNA from PCR and other enzymatic reactions
Technology	Silica technology
Type(s) of DNA recovered	ss DNA and dsDNA

<http://www.qiagen.com/products/catalog/sample-technologies/dna-sample-technologies/dna-cleanup/qiaquick-pcr-purification-kit#technicalspecification>

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Ligation Protocol with T4 DNA Ligase (M0202)

Protocol

- Sat up the following reaction in a microcentrifuge tube on ice.
(T4 DNA Ligase should be added last. Note that the table shows a ligation using a molar ratio of 1:3 vector to insert.)

COMPONENT	20 µl REACTION
10X T4 DNA Ligase Buffer*	2 µl
Vector DNA (3 kb)	50 ng (0.025 pmol)
Insert DNA (1 kb)	50 ng (0.076 pmol)
Nuclease-free water	to 20 µl
T4 DNA Ligase	1 µl

* The T4 DNA Ligase Buffer should be thawed and resuspended at room temperature.

- Gently mix the reaction by pipetting up and down and microfuge briefly.
- For cohesive (sticky) ends, incubate at 16°C overnight or room temperature for 10 minutes.
- For blunt ends or single base overhangs, incubate at 16°C overnight or room temperature for 2 hours (alternatively, high concentration T4 DNA Ligase can be used in a 10 minute ligation).
- Chill on ice and transform 1-5 µl of the reaction into 50 µl competent cells.

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